



# Undulator Spectra as a Diagnostic for Electron Beam Energy Spread and Emittance

M. Bakeman, J. van Tilborg, T. Sokollik,  
C. Toth, W. Fawley, C. Schroeder, K. E.  
Robinson, E. Esarey, and W. Leemans

Lawrence Berkeley National  
Laboratory



# Introduction



---

- Loasis 1GeV Laser Wakefield Accelerator
- Undulator Radiation as a beam diagnostic
- THUNDER Undulator Characterization
- Quadrupole Electron Beam Transport
- XUV Spectrometer Design and Calibration
- Conclusions and Future Work

- Laser:  $a_0 \sim 1.46$  (40 TW, 37 fs)
- Capillary:  $D = 312 \mu\text{m}$ ;  $L = 33 \text{ mm}$



1 GeV beam




	Sim	Expt
Q (pC)	25-60	35
E (GeV)	1.0	1.1
dE/E RMS (%)	4	2.5
div. (mrad)	2.4	1.6

# Undulator Radiation Diagnostic to Measure Emittance and Energy Spread


- Single shot electron beam measurements:
  - charge, energy, energy spread, transverse emittance

## Schematic layout:





- Supporting systems:
  - Quads for refocusing
  - BPMs for positioning
  - OTR, phosphor screens

# Emittance and energy spread affect the undulator spectra



Energy spread  $\times 10$



Emittance  $\times 10$







# Undulator radiation based diagnostic will allow for single shot $\Delta E/E$ and $\varepsilon$ measurement


## Changing energy spread



## Line Width vs. Energy Spread



# Experimental Setup in the Loasis lab





## Magnetic Characteristics



- Length = 1.5m
- Magnetic period = 2.18cm
- Number of periods = 66
- Peak Field = 0.61T
- Wiggler Parameter  $K = 1.25$
- Betatron period = 3.7m









# THUNDER tuned and characterized with the Spring-8 SPECTRA code




E-beam trajectory before tuning




E-beam trajectory after tuning



# Quadrupoles transport E-beam to undulator




# New capillary chamber designed to fit quads and stages



# Micro-Channel Plate Based Spectrometer from 8 - 62nm


## Optical Characteristics

- First harmonic = 20nm
- Natural divergence = 1 mrad
- Pulse length = < 25 fs
- First harmonic flux =  $1 \times 10^5$




## XUV Spectrometer properties

- 8-62nm spectral range
- Four jaw Aperture System
- MCP detector with CsI photocathode
- Aberration corrected concave grating






# Characterization completed at the Advanced Light Source




- Characterized on beamline 6.3.2
- 8-62nm Spectral Range
- 0.07nm resolution at 20nm
- Hitachi Aberration Corrected Concave Grating




# XUV Spectrometer Calibration

## Spectral Resolution vs. Wavelength



## Wavelength vs. Detector Position







# Conclusions and Future Work

---

## First experiments to begin this spring at Loasis

- Undulator radiation can be used as a single shot emittance and energy spread beam diagnostic
- Experiments will be conducted using the Loasis LWFA at 500MeV
- Undulator characterized and installed on beamline
- XUV spectrometer built and characterized
- Quadrupoles and magnetic spectrometer will be installed this spring

This work is supported by the Defense Threat Reduction Agency and the Director, Office of Science, High Energy Physics, U.S. Dept. of Energy under Contract no. DE-AC02-05CH11231